



SUN RISE
EYE CARE PVT. LTD.

3 & 4, Karan Markt, 1st Floor, Opp. Gopal Mandir,
Panj Peer, Jalandhar-144 001.
Tel. : 2211837 Fax : 91-181-2211837
E-mail : srecpl@jla.vsnl.net.in

SRECPL/4685/2014-015.
1st October,2014.

Rotary Club, North Ludhiana,
A/C Viaayanand Diagnostic Centre,
Basant Road, Civil Lines,
LUDHIANA-141001.

Dear Sir,

As per telephonic discussion with the under of Dr. Sameer Dogra Medical Suptd. of Viayanand Diagnostic Centre, regarding the Ophthalmic Equipments of TOPCON AND PLANTECH and thanks for interest shown in TOPCON AND PLANTECH products.

As desired by you we are pleased to enclosed our Quotations with Literature for the items required by you for your kind perusal.

In case if you require any further information/clarification or any assistance please let us know and the undersigned will call on you personally.

Thanking you.

Yours faithfully,
For SUN RISE EYE CARE PVT. LTD.


RAVINDER LUTHER.

Encl:As above.

Quotation No. **3014**

Dated 1.10.2014.

ROTARY CLUB NORTH. LUDHIANA.
A/C, VIJAYANAND DIAGNOSTIC CENTRE,
BASANT ROAD, CIVIL LINES,
LUDHIANA.141001.



Sun Rise

EYE CARE PVT. LTD.

3 & 4, Karan Market, 1st Floor, Opp. Gopal Mandir,
Panj Peer, Jalandhar City - 144 001

Ph.: 2211837, 5000837, Fax.: 91-181-2211837, Mob.: 98144-40837

E-mail : sreclp@jla.vsnl.net.in

CURRENCY OF SALE

INDIAN RUPEES


QNTY.	DESCRIPTION	RATE	PER	AMOUNT
	<u>TOPCON OPHTHALMIC INSTRUMENTS.</u>			
ONE NOS.	TOPCON SLIT LAMP SL-2G.	2,10,000.00	SET	<u>2,10,000.00</u>
ONE NOS.	PLANTECH MOTORISED INST.TABLE.	15,000.00		<u>15,000.00</u>

For **SUN RISE EYE CARE PVT. LTD.**



Managing Director / Director

 **TOPCON**
CONNECTING VISIONS

TERMS AND CONDITION APPLICABLE
FOR THIS PROFORMA

1. PAYMENT : 100% Advance by demand draft in MEHRA EYETECH PVT. LTD.,
MUMBAI, along with confirmation of order.
Prices inclusive of PONDICHERY PGST & CST
are subject to PONDICHERY Jurisdiction.
2. DELIVERY : 4 weeks after receipt of order with advance.
3. PRICES : Prices quoted are CIF LUDHIANA all other charges like freight,
Insurance, packing etc. are inclusive.
4. VALIDITY :  30 Days.
5. COUNTRY OF ORIGIN : JAPAN MIT INDIA.
6. INSTALLATION : Installation charges are inclusive.
7. UNPRECEDEN TED EVENTS : All prices quoted will be subject to change without notice in event of
unprecedented events like war, fresh government measures pertaining
to imports, price would be applicable as per the prevalent price at the
time of delivery.
8. GUARANTEE : Equipments are under warrantee for a period of one year from the date
of supply.
9. DISPUTES : Any disputes will be subject to PONDICHERY jurisdiction.
10. Force Majure Clause : Quotations are agreements to supply goods as per the proforma is
subject to force majure clause.
11. Octroi / Entry Tax : If applicable will be extra.

Slit Lamp SL-2G



Basic, Economical, Environmentally Friendly

Long Life LED Illumination Saves Power, Helps the Environment and Reduces Bulb Replacement.



Features


- » Long Life and Power-Saving LED Illumination
- » Easy to Install with Incorporated AC Power Supply
- » Three Magnifications (10x; 16x; 25x)
- » Compact and Practical , Ergonomic Design

- SL-2G LED Slit Lamp

Still Using Halogen Lamp Illumination?


Light Emission Diodes (LEDs) are increasingly being applied to optical and medical instrumentation thanks to their clear and bright light, low heat illumination and long life. LEDs also have the added advantage of being environmentally friendly reducing the waste of frequently replacing halogen and tungsten bulbs.

As a leader in Ophthalmic Instrumentation, Topcon is incorporating LED as the main illumination in its SL-2G slit lamp. The approx 10,000 hours lifetime of the LED is expected of 100 times longer than conventional halogen bulbs. This makes the SL-2G easier to maintain, economical and eco-friendly.



Quality of Illumination

Uniform and consistent illumination is achieved by LED source.



Three Magnifications

10x, 16x, 25x cover the needs for a comprehensive slit lamp examination.



Smooth Illumination Adjustment

The LED illumination can be adjusted continuously from 0 to maximum.



AC Adapter

The AC adapter included with the SL-2G enables the use of any power source and adds flexibility to the installation in different types of furniture and delivery systems.



Specifications

Microscope

Type	Galileo convergence
Magnification selection	3 steps by rotation drum
Magnification	10x, 16x, 25x
Eyepiece	12.5x
Diopter adjustment	-5D ~ +3D
Total magnification (field of view)	10x (φ 22.5 mm) 15.98x (φ 14.1 mm) 25.53x (φ 8.8 mm)

Illumination

Slit width	Continuous from 0 to 14 mm
Slit length	Continuous from 1 to 14 mm
Aperture	Continuous from 1 to 14 mm , φ 0.3 , 1 , 5 , 10 , 14 mm
Slit angle	0° to 180°
Filter	Blue, Redfree
Light source	LED

Power

Power consumption	40VA
Input (Primary)	AC 100-240V 50-60Hz

Dimensions and Weight

Standard (with table top and chinrest)	550mm(W) x 399mm(D) x 520mm(H), 15.2 Kg
----------------------------------------	-----------------------------------------

Optional Accessories



» External Fixation Target (SO-FT02)



» Hruby Lens (SO-HL01)



» Yellow Filter



» Tonometer Mount for R900 (SO-TM)



» Tonometer Mount for 870 (SO-TM2)

- » 12.5x Measuring Eyepiece
- » 20X Eyepiece
- » Parallel Binocular Tubes (PB-2)



SUN RISE
EYE CARE PVT. LTD.

3 & 4, Karan Market, 1st Floor, Opp. Gopal Mandir,
Panj Peer, Jalandhar-144 001
Tel. : 2211837, 5000837 Tele/Fax : 91-181-2211837
Mob. : 98144-40837
E-mail : srepl@jia.vsnl.net.in

*Subject to change in design and/or specifications without advanced notice.

IMPORTANT

In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.



TOPCON CORPORATION



TOPCON CORPORATION



TOPCON CORPORATION



TOPCON CORPORATION

75-1 Hasunuma cho, Itabashi-ku, Tokyo 174-8580, Japan Phone 3.3558-2523/2522 Fax 3.3965-6898 www.topcon.co.jp

TOPCON MEDICAL SYSTEMS, INC.

111 Bauer Drive Oakland, NJ 07436, U.S.A.
Phone +1-201-599-5100 Fax +1-201-599-5250 www.topconmedical.com

TOPCON MEDICAL LASER SYSTEMS, INC.

3130 Coronado Drive Santa Clara, California 95054 USA
Phone +1-408-235-8200 Phone (USA only) +1-888-760-8657
Fax +1-408-235-8259 E-mail: mlmsinfo@topcon.com

TOPCON CANADA INC.

110 Provencer Avenue, Boisbriand, QC J7G 1N1 CANADA
Phone +1-450-430-7771 Fax +1-450-430-6457 www.topcon.ca

TOPCON EUROPE MEDICAL B.V.

Esplanade 11, 2908 LJ Capelle a/d IJssel, P.O. Box 145,
2900 AC Capelle a/d IJssel, THE NETHERLANDS
Phone +31-(0)10-4585077 Fax +31-(0)10-4585045
E-mail: medical@topcon.nl www.topcon.eu

ITALY OFFICE

Viale dell'Industria 60, 20037 Paderno Dugnano, (Milano), ITALY
Phone +39-02-9186671 Fax +39-02-91081091 E-mail: topconitaly@iscall.it www.topcon.it

DANMARK OFFICE

Præstemarkvej 25, 4000 Roskilde, DANMARK
Phone +45-46-327500 Fax +45-46-327555
E-mail: topcon@topconandmark.dk www.topconandmark.dk

IRELAND OFFICE

Unit 276, Blanchardstown Corporate Park 2 Ballycoolin Dublin 15, IRELAND
Phone +353-18975900 Fax +353-18293915 E-mail: medical@topcon.ie www.topcon.ie

TOPCON S.A.R.L.

BAT A1 3 route de la révoite 93206 SAINT DENIS CEDEX
Tel +33 1 49 21 23 23 Fax +33 1 49 21 23 24 E-mail: topcon@topcon.fr www.topcon.fr

TOPCON DEUTSCHLAND G.m.b.H.

Hanns-Martin-Schleyer Strasse 41, D-47877 Willich, GERMANY
Phone +49-(0)2154-8850 Fax +49-(0)2154-885177 E-mail: med@topcon.de www.topcon.de

TOPCON SCANDINAVIA A.B.

Noongalan 2, P.O. Box 25 43151 Mölndal, SWEDEN
Phone +46-(0)31-7199200 Fax +46-(0)31-7109249 E-mail: medical@topcon.se www.topcon.se

TOPCON ESPAÑA S.A.

HEAD OFFICE
Frederic Mompou 4 Esc. A Bajos 3, 08960 Sant Just Desvern Barcelona, SPAIN
Phone +34 93-4734057 Fax +34-93-4733932 E-mail: medical@topcon.es www.topcon.es

PORTUGAL OFFICE
Rua da Forte, 6-6A L. O 22 2790-072 Carnaxide, PORTUGAL
Phone +351 210-994626 Fax +351 210-938786 www.topcon.pt

TOPCON (GREAT BRITAIN) LTD.

Topcon House, Kennet Side Bono Lane Newbury Berkshire RG14 5PX United Kingdom
Phone +44-(0)1635-551120 Fax +44-(0)1635-551170 E-mail: medical@topcon.co.uk www.topcon.co.uk

TOPCON POLSKA Sp. z o.o.

ul. Warszawska 23, 42-470 Siewierz, POLAND
Phone +48-(0)32-6705045 Fax +48-(0)32-6713405 www.topcon.polska.pl

TOPCON SINGAPORE MEDICAL PTE. LTD.

60 Alexandra Terrace, #08-27 The Comtech, Singapore 118502
Phone +65-68720606 Fax +65-67736150 E-mail: medical_sales@topcon.com.sg www.topcon.com.sg

TOPCON INSTRUMENTS (MALAYSIA) SDN.BHD.

No D1 (Ground Floor), Jalan Exzodia 2, Off Jalan Ampang Putra,
Taman Ampang Hilir, 55100 Kuala Lumpur, MALAYSIA
Phone +60-(0)3-42709866 Fax +60-(0)3-42709766

TOPCON INSTRUMENTS (THAILAND) CO.,LTD.

77/162 Sinsathorn Tower, 37th Floor, Krungthornburi Rd. Klongtonnai,
Klongsarn, Bangkok 10500, THAILAND
Phone +66(0)2-440-1152-7 Fax +66(0)2-440-1158

TOPCON CORPORATION BEIJING OFFICE

Block No 9, Kangding Street Beijing Economic Technological Development Area,
Beijing 100176, CHINA
Phone +86-(0)10-6780-2799 Fax +86-(0)10-6780-2790

TOPCON CORPORATION SHANGHAI OFFICE

14L Huamin Empire Plaza, No 726, Yan'an Xi Road,
Shanghai 200050, CHINA
Phone +86-(0)21-5238-7722 Fax +86-(0)21-5237-0761

TOPCON CORPORATION BEIRUT OFFICE

P.O. Box 70-1002 Antlias Borut, LEBANON
Phone +961-4-523825/523526 Fax +961-4-521119

TOPCON CORPORATION DUBAI OFFICE

P.O. Box 293705 Dubai Airport Free Zone L U J-12, Dubai, U.A.E.
Phone +971-4-293-5900 Fax +971-4-299-5901