

QUOTATION OFFER FOR

INFANT INCUBATOR OF DELIVERY, INSTALLATION & BASIC USER TRAINING

TO
THE PRESIDENT-ELECT
ROTARY CLUB OF PORT-HARCOURT
DISTRICT 9140
RIVER STATE

MEDICAL EQUIPMENT



Scooda Enterprises Limited

Mobile: +2348060088969

Email: Scoodaenterprisesltd@gmail.com and dedeoghon@gmail.com

......Designed for your Experience !!!



Client: The President-Elect

Rotary Club Of Port-Harcourt District 9140

River State.

Ref No: SE-02112016

Currency: Naira

Date: 2nd November 2016

Manufacturer: Bistos, Germany & South Korea.

Dear Sir

We are pleased to present our offer for Medical Equipment / Consumables for your kind consideration.

S/N	PICTURE	NAME OF ITEM	Manfacturer	Qty	UNIT PRICE	TOTAL PRICE
	INFANT INCUBATOR					
1	BISTOS BT-500 INFANT INCUBATOR	Cotroles Advantage A	BISTOS BT-500 is an infant incubator, which can precisely control temperature and humidity with low noise providing the best environment for the recovery of neonate. Also, CCD Camera enables remote monitoring of neonate. Also, CCD Camera enables remote monitoring of neonate. Based on ergonomic design. BT-500 will provide extreme comport to customers. FEATURES • Stable / air skin temperature control • Accurate humidity control • Comfortable & soft tilting structure • 7" color TFT LCD display • Powerful lifting structure (optional) • OZ monitoring (optional) • Masimo SpO2 technology (optional) • Masimo SpO2 technology (optional) • Physical characteristics (with standard accessories) • Size 140 x 103 x 141 cm (H x L x D) • Weight 89 kg • Temperature Control • Air temperature control range 23.0 ~ 37.0 ± 0.3 °C (Override >39.0 °C) • Skin temperature control range 35.0 ~ 37.5 ± 0.3 °C (Override >39.0 °C) • Peripheral temperature YES (Option) • Humidity Control • Humidity control range 40 ~ 95% ± 5 % RH • Measurement rage 15 ~ 99% ± 5 % RH • Measurement rage 15 ~ 99% ± 5 % RH • Measurement rage 15 ~ 99% ± 5 % RH • Measurement rage 15 ~ 99% ± 5 % RH • Measurement page 17" TFT COLOR LCD • Trend up to 7 days • Alarm 19 kind of alram • Multi-language support YES • Hood • Hood size 39 x 91 x 51 cm (H x L x D) • Matress size 38 x 73 cm • Matress tilt 12 • Air Velocity < 10 cm/s (at 10cm above the center of mattress) • Warming-up time (from ambient 22 °C) < 35 min. • Noise level < 45 dBA • Micro air filter perficiency 99.8 % • Option • Fixed Stand(FX) • Dimension 100cm(W) x 81cm(D) x 80cm(H) • Basket 41cm x 42cm, 10kg(Weight Limit) • Lifting Stand(LS) • Dimension 100cm(W) x 81cm(D) x 62cm~82cm(H) • Speed 180mm/min • Basket 41cm x 42cm, 20kg(Weight Limit) • Weighting Scale • Measurement 0~10kgf, 50g • Masimo SpO2 • SpO2 1~100% • Pulse rate 25bpm ~ 240bpm • Accuracy - SpO2: Nenate 70~100% ± 3 digits • Pulse rate 25 c ~ 240 ± 3 digits • Monitoring System External 7" color TFT LCD • CCD Camera :510 x 492 (Resolution	1	4,863,936.00	4,863,936.00
		•			Sub-Total	4,863,936.00

Sub-Total	4,863,936.00			
Delivery,Installation and Basic UserTraining	250,000.00			
Total Value of the Offer:	5,113,936.00			



Invoice

Client: The President-Elect Rotary Club Of Port-Harcourt

District 9140 River State. Ref No: SE-02112016

Currency: Dollars Ex: N 470

Date: 2nd November 2016

Manufacturer: Bistos, Germany & South Korea.

Dear Sir

We are pleased to present our offer for Medical Equipment / Consumables for your kind consideration.

	We are pleased to present our offer for Medical Equipment / Consumables for your kind consideration.					
S/N	PICTURE	NAME OF ITEM	Manfacturer	Qty	UNIT PRICE	TOTAL PRICE
	INFANT INCUBATOR					
1	BISTOS BT-500 INFANT INCUBATOR	Dodrvich Management of the state of the stat	BISTOS BT-500 is an infant incubator, which can precisely control temperature and humidity with low noise providing the best environment for the recovery of neonate. Also, CCD Camera enables remote monitoring of neonate, and in-built SpO2 monitor and weighing scale provide total monitoring of neonate. Based on ergonomic design. BT-500 will provide extreme comport to customers. FEATURES Stable / air skin temperature control Accurate humidity control Comfortable & soft tilting structure "" color TFT LCD display Powerful lifting structure (optional) O2 monitoring (optional) Masimo SpO2 technology (optional) SPECIFICATION Physical characteristics (with standard accessories) Size 140 x 103 x 141 cm (H x L x D) Weight 89 kg Temperature Control Air temperature control range 23.0 ~ 37.0 ± 0.3 °C (Override >39.0 °C) Skin temperature control range 35.0 ~ 37.5 ± 0.3 °C (Override >39.0 °C) Peripheral temperature YES (Option) Humidity Control Humidity Control range 40 ~ 95% ± 5 % RH Measurement rage 15 ~ 99% ± 5 % RH Control system Ultrasonic & steam Water tank Capacity 1,000 ml Display Display panel 7" TFT COLOR LCD Trend up to 7 days Alarm 19 kind of airam Multi-language support YES Hood Hood size 39 x 91 x 51 cm (H x L x D) Matress tilt 12 Air Velocity < 10 cm/s (at 10cm above the center of mattress) Warming-up time (from ambient 22 °C) < 35 min. Noise level < 45 dBA Micro air filter particle size 0.3 micron Micro air filter particle size 0.3 micron Micro air filter particle size 0.3 micron Micro air filter afficiency 99.8 % Option Fixed Stand(FX) Dimension 100cm(W) x 81cm(D) x 80cm(H) Basket 41cm x 42cm, 10kg(Weight Limit) Lifting Stand(LS) Dimension 100cm(W) x 81cm(D) x 62cm~82cm(H) Speed 180mm/min Basket 41cm x 42cm, 10kg(Weight Limit) Weighting Scale Measurement 0~10kgf, 50g Masimo SpO2 SpO2 1~100% Pulse rate 25 pbm ~ 240bpm Accuracy -SpO2: Nenate 70~100% ± 3 digits Monitoring System Extend 7" Color TFT LCD CCD Camera :510 x 492 (Resolution) Masimo techn	1	10,348.80	10,348.80
					Sub-Total	10 3/8 80

Sub-Total	10,348.80
Delivery,Installation and Basic UserTraining	500.00
Total Value of the Offer:	10,848.80



TERMS & CONDITIONS:

Ref: #REF!

Dear Madam,

We are pleased to submit to you our offer according to your specification. We are enclosing our proforma invoice No:SE-02112016 of 2.11.2016.

Country of Origin: Japan

Delivery Time: A complete delivery approx. within 2 - 3 weeks after receipt of order and payment.

Mode of payment: 70% Advance payment of the total sum,30% of balance paid after Delivery &

Installation.

Spare Parts: Spare parts availability for ten years from date of installation.

Warranty: Our standard manufacturer's warranty is 12 months from date of

Please note, that in the event of any variation to the above terms and conditions, or quantity of products as quoted, we reserve the right to review our offer.

We hope our offer meets your interest, and look forward to hearing from you.should you have any further questions, please do not hesitate to contact us.

Best regards

Davidson Edeoghon

Head Marketing and Strategy

Scooda Enterprises Limited Mobile: +2348060088969

Email: Scoodaenterprisesltd@gmail.com or dedeoghon@gmail.com

......Designed for your Experience !!!